

**GLUT-1 (SPM498)**

Catalog # MS-10637-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200ug/ml)

Catalog # MS-10637-R7 (7.0ml)

Catalog # MS-10637-PCS

**INTENDED USE:**

- **For In Vitro Diagnostic Use:** This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.
- **Description:** Glucose is fundamental to the metabolism in mammalian cells. Several glucose transporter protein (Glut) isoforms have been identified and shown to function in response to insulin and IGF-1 induced signaling. GLUT-1 is detectable in many human tissues including those of the colon, lung, stomach, esophagus, and breast. GLUT-1 immunoreactivity in some cancers, including trans carcinoma of the urinary bladder, has been associated with aggressive behavior.
- **Expected Staining Pattern:** Cell membrane
- **Positive Control:** Esophageal CA

**MATERIALS PROVIDED:**

**GLUT-1 (SPM498) (refer to catalog number):**

- #MS-10637-P (or -P0, -P1) 200ug/ml of antibody purified from ascites. Prepared in 10mM PBS, pH 7.4, with 0.2% BSA and 0.09% sodium azide.
- or
- #MS-10637-R7: (7.0ml) of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide.
- or
- #MS-10637-PCS: 5 positive control slides.
- **Antibody Concentration:** 200ug/ml
- **Host:** Mouse
- **Epitope:** c-terminal
- **Species Reactivity:** Human and Rat. Others not tested.
- **Clone Designation:** SPM498
- **Ig Isotype / Light Chain:** IgG1/k
- **Immunogen:** A synthetic peptide derived from C-terminal of human GLUT-1
- **Microbiological State:** This product is not sterile.

**MATERIALS REQUIRED, BUT NOT PROVIDED:**

- **Antibody Diluent:** For concentrated antibodies, the antibody must be diluted before using. Use Lab Vision Antibody Diluent (catalog # TA-125-UD). Refer to diluent product instructions for use.
- **Negative Control Reagent:** Refer to the "General Protocol" instructions.
- **Visualization System:** Refer to the "General Protocol" instructions.

**METHODS AND PROCEDURES:**

<b>Specimen Preparation</b>	Refer to the "General Protocol" instructions.
<b>Dilution of Concentrated Antibody</b>	1:200 in antibody diluent
<b>Tissue Section Pretreatment</b>	Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM citrate buffer, pH 6.0 (Lab Vision catalog # AP-9003), for 10-20 minutes followed by cooling at room temperature for 20 min..
<b>Primary Antibody Incubation Time</b>	30 mins at Room Temperature
<b>Visualization</b>	To detect antibody, follow the instructions provided with the visualization system.

**STORAGE and STABILITY:**

This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.