

ALK (Anaplastic Lymphoma Kinase) / p80 Ab-1

Catalog # MS-1104-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)

Catalog # MS-1104-R7 (7.0ml)

Catalog # MS-1104-PCS

INTENDED USE:

- **For In Vitro Diagnostic Use:** This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.
- **Description:** ALK / p80 is a hybrid of the anaplastic lymphoma kinase (ALK) gene and the nucleophosmin (NPM) gene resulting from the t(2;5)(p23;q35) translocation found in a third of large cell lymphomas.
- **Expected Staining Pattern:** Cytoplasmic and/or nuclear
- **Positive Control:** Anaplastic lymphoma

MATERIALS PROVIDED:

ALK (Anaplastic Lymphoma Kinase) / p80 Ab-1 (refer to catalog number):

- #MS-1104-S (or -S0, -S1): Tissue culture supernatant, concentrated, with 0.09% Sodium Azide.
or
- #MS-1104-R7: (7.0ml) of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide.
or
- #MS-1104-PCS: 5 positive control slides.
- **Antibody Concentration:** Not known
- **Host:** Mouse
- **Mol. Wt. of Antigen:** 80kDa
- **Epitope:** aa 419-520 (tyrosine kinase domain)
- **Species Reactivity:** Human. Others not-tested.
- **Clone Designation:** 5A4
- **Ig Isotype / Light Chain:** IgG1
- **Immunogen:** Recombinant protein corresponding to a region which spans the tyrosine kinase catalytic domain and part of the C-terminus of the NPM-ALK transcript (419-520aa).
- **Microbiological State:** This product is not sterile.

MATERIALS REQUIRED, BUT NOT PROVIDED:

- **Antibody Diluent:** For concentrated antibodies, the antibody must be diluted before using. Use Lab Vision Antibody Diluent (catalog # TA-125-UD). Refer to diluent product instructions for use.
- **Negative Control Reagent:** Refer to the "General Protocol" instructions.
- **Visualization System:** Refer to the "General Protocol" instructions.

METHODS AND PROCEDURES:

Specimen Preparation	Refer to the "General Protocol" instructions.
Dilution of Concentrated Antibody	1:40-1:80 in antibody diluent
Tissue Section Pretreatment	Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM citrate buffer, pH 6.0 (Lab Vision catalog # AP-9003), for 10-20 minutes followed by cooling at room temperature for 20 min.
Primary Antibody Incubation Time	60 minutes at Room Temperature
Visualization	To detect antibody, follow the instructions provided with the visualization system.

STORAGE and STABILITY:

This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.

REFERENCES:

- 1) Pittaluga S, et al. American Journal of Pathology. 151(2):343-351 (1997).
- 2) Pulford K, et al. Blood. 89(4): 1394-1404 (1997)
- 3) Downing J R, et al. Blood. 85(12): 3416-3422 (1995).