

Large Volume UltrAb Diluent Plus

INTENDED USE

For In Vitro Diagnostic Use

<u>AVAILABILITY:</u>	<u>Catalog #</u>	<u>Volume</u>
	TA-125-UDX	125 ml

<u>SPECIFICITY:</u>	N/A
<u>ENZYME:</u>	N/A
<u>CHROMOGEN/SUBSTRATE:</u>	N/A

DESCRIPTION

UltrAb Diluent is prepared for use in immunochemical procedures. It is especially useful with mouse IgG1 and rabbit antibodies, resulting in stronger staining. Antibodies can be diluted directly in this reagent and can be stored for up to 18 months when refrigerated (extremely high dilutions will result in reduced shelf life). Dilute primary antibodies, or biotinylated antibodies in this reagent.

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at 2-8°C. Product is stable for 24 months from the date of manufacture.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

Dilute antibody directly in this reagent.

REFERENCES

N/A

TROUBLESHOOTING

Please contact Thermo Fisher Scientific Technical Support by phone (1-510-991-2800 or 1-800-828-1628) or by email (lab.reagents@thermofisher.com).